Soffio Cube

INSTRUCTION MANUAL • MODE D'EMPLOI • INSTRUCCIONES DE USO • INSTRUÇÕES DE USO
GEBRUIKSAANWIJZING • INSTRUKCJA OBSŁUGI • BRUKSANVISNING • NÁVOD K POUŽITÍ
ΟΔΗΓΙΕΣ ΧΡΗΣΕΩΣ • مَادِخْتَسَالَا تَامِ لِعَت

 **Air Liquide**
HEALTHCARE

ENGLISH • FRANÇAIS • ESPAÑOL • PORTUGUÊS • NEDERLANDS • POLSKI •
SVENSKA • ČEŠTINA • ΕΛΛΗΝΙΚΑ • دِي بَرَعَالَا

Fig. A

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
Fig. B

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

Fig. C1

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

Fig. C2

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Fig. D

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1.3 - IMPORTANT NOTES

- Read the instructions carefully and keep them in a safe place for future reference.
- Do not use the device or accessories for anything other than the intended purpose.
- Children must only be allowed to use the device under the supervision of an adult who has read and understands the instructions.
- Some components are small enough to swallow and may lead to suffocation.
- The power cable and connecting pipe are long enough to pose a strangulation hazard.
- Read the list of components in the table in section 7 to see if there are any components that have caused an allergic reaction in the past.
- During the session, the device must stand on a firm flat surface.
- Never cover the air vents during operation.
- Keep well away from sources of heat.
- Do not touch the plug or the electrical components wet or with damp hands.
- Device protected against finger access to hazardous parts and protected against vertically falling drops (IP21).
- Do not pour water or other liquids onto the device. Prevent water and other liquids from getting in through the air vents.
- Do not use in the bath or shower.
- If the device should fall into liquid, unplug before removing it. Do not use the device. Send it to an approved service centre for an overhaul.
- Beware of compressed air coming out of the air outlet as this may cause a hazard; it must only be used to supply the nebuliser.
- The device is not suitable for use in the presence of anaesthetic mixtures with air, oxygen or nitrogen protoxide.
- Performance of the device is only guaranteed if all the accessories used are original ones. Only original spare parts must be used.
- Do not use the device if the electrical cable is damaged.
- Check that the casing is integral and there is no sign of dents or damage.
- Never open the device.
- Never bend the nebuliser over 60°.
- The manufacturer declares that the device has a useful life of at least 5 years.

1.4 - DESCRIPTION OF THE DEVICE

The SOFFIO CUBE aerosol therapy device is made up of a compressor driven by an electric motor in the basic unit and a complete accessory kit.

1.5 - PATIENT GROUPS

The SOFFIO CUBE aerosol therapy device provides highly effective inhalation therapy for patients of any age, from young children to the elderly.

1.6 - CONTENTS

The package contains (Fig. A and C):

- The basic unit (1) with power cable (2), ON/OFF switch (3), air filter slot and filter (4), air intake (5) and nebulizer support (6).
- An accessory kit composed by: nebuliser (upper body (7), internal atomiser (8), lower body (9)), mouthpiece (10), nosepiece (11), adult mask (12), pediatric mask (13), connecting tube (14), connector (15)
- Replacing air filters
- Instructions for use
- Bag

2 - PREPARATION AND USE

Before use, make sure the voltage indicated on the rating plate underneath the device corresponds to the mains voltage and the wiring system complies with the current regulations.

Prior to using the device for the first time, we recommend cleaning it as described in the section «Cleaning, disinfecting and sterilizing».

1. Fill the nebuliser with the inhalation solution as per your doctor's instructions (Fig.B). Ensure that you do not exceed the maximum level.
2. Assemble the nebuliser and connect the mouthpiece or the nosepiece or the mask. (Fig.C). Ensure that all parts are complete.
3. Connect the nebuliser with the connecting tube (14) to the compressor (1) and plug the power cable (2) into the socket.
4. To turn on the device, push the ON/OFF switch (3) onto the position «I». Lean the mouthpiece between the lips or the nosepiece on the nostrils or lean the mask to the face over mouth and nose.
5. Sit in a relaxed position with the upper body upright. Breathe in and out calmly during the therapy. Do not lie down while inhaling. Stop inhalation if you feel unwell.
6. After completing the inhalation period recommended by your doctor, switch the ON/OFF switch (3) to position «O» to turn off the device and unplug it from the socket.
7. Empty the remaining medication from the nebuliser and clean the device as described in the section «Cleaning, disinfecting and sterilizing».

3 - CLEANING, DISINFECTING AND STERILIZING

Always wash your hands well before cleaning and disinfecting the accessories.

Do not expose the compressor to water or heat.

Replace the connecting tube for each treatment with a new patient or in case of impurities.

Unplug the device before cleaning.

Cleaning

The compressor unit (1) and the connecting tube (14) should be cleaned with a clean, moist cloth.

Clean all the nebuliser components and its accessories under warm tap water (max. 60°C) adding if necessary a small quantity of detergent following dosage and use limitations as provided by detergent manufacturer. Rinse thoroughly making sure that

all deposits are washed away and leave to dry.

Disinfecting

The compressor unit (1) can be disinfected with a clean cloth and some chemical disinfectant following dosage and use limitations as provided by disinfectant manufacturer.

All the accessories (except the connecting tube) can be disinfected with chemical disinfectants following dosage and use limitations as provided by disinfectant manufacturer. Disinfectants are usually available at pharmacies.

Steam sterilizing

The nebuliser and its accessories (except air hose and masks) can be heat steam sterilized up to 121 °C (20 min.) or 134 °C (7 min.). EN554/ISO11134. The sterilization packaging must conform to EN868/ISO11607 and be suitable for steam sterilization. After sterilization always let all components cool down to ambient temperature before further use.

Do not repeat sterilization cycle when components are still warm.

4 - MAINTENANCE

The SOFFIO CUBE aerosol therapy device does not require any particular maintenance. The compressor does not need to be lubricated as the piston is the dry-operating type. It is advisable replace nebulizer after some 100 to 120 treatments on single patient or after about 20 sterilization cycles. Check the filter continually for cleanliness and replace it if dirty, or after a maximum of 3 months use. Spare filters are provided with the device. To replace the filter, open the air filter slot (4) (if needed using a tool) and replace the filter (Fig. D).

5 - PROBLEMS, POSSIBLE CAUSES AND SOLUTIONS

The device cannot be switched on	Ensure the power cord (2) is correctly plugged into the socket. Ensure the ON/OFF switch (3) is in the position «I».
The nebuliser functions poorly or not at all	Ensure that the connecting tube (14) is correctly connected at both ends. Ensure the connecting tube is not squashed, bent, dirty or blocked. If necessary, replace with a new one. Ensure the nebuliser is correctly assembled. Ensure the required medication has been added and in the appropriate quantity. Ensure that the small hole of the nebulizer's lower body (9) is not clogged.

If you have failed to solve the problem after performing the above checks and the device still does not work, have it checked by an authorized service centre.


6 - TECHNICAL DATA

Supply voltage	230 V ~ 50 Hz
Operation	continuous

Maximum air pressure at the device	1,85 bar
Maximum air flow at the device	10 l/min
Operating air flow with the nebulizer	5 l/min
Nebulizing rate (with NaCl 0,9% solution)	0,35 ml/min
Maximum absorbed power	150 VA
Dimensions of basic unit	135x140x155H mm
Weight of the basic unit	1,1 kg approx
Noise level at 1 m	58 dBA
Minimum nebulizer filling volume	2ml
Maximum nebulizer filling volume	12 ml
Aerosol output*	81,2 µl
Aerosol output rate*	30,2 µl/min
MMAD* (Graph 1)	2,94 µm

* measures performed as indicated in UNI EN 13544-1:2009

Graph 1

The y-axis of the graphs represents the percentage of volume of atomised solution with sizes lower than the corresponding dimensional value expressed in µm on the x-axis.

50% by volume of the atomised particles is less than the MMAD (Median Mass Aerodynamic Diameter). This value is indicative of the Nebuliser's ability to produce particles suitable for a specific therapeutic treatment.

The shown MMAD value has been measured with a cascade impactor NGI Next Generation Impactor.

7 - DISPOSING OF THE DEVICE AND ITS COMPONENTS

At the end of their useful life, the device and its components must be disposed of in accordance with the current regulations and environmental laws. If no legal obligations exist, the various materials must be sorted and disposed of separately. The table below shows the materials of which the various components are made.

None of the materials used contains phthalates or traces of natural latex.

COMPONENT	MATERIAL
Compressor unit	Waste from Electric and Electronic Equipment (WEEE)
External casing of the compressor unit	Polypropylene
Nebulizer and connector	Polypropylene
Mouthpiece	Polypropylene
Nosepiece	Polypropylene
Masks	Non toxic PVC
Connecting tube	Non toxic PVC
Box	Cardboard
Instructions for use	Paper
Bag	Nylon
Air filters	Needled polypropylene

8 - ELECTROMAGNETIC COMPATIBILITY

This device complies with the requirements of EN 60601-1-2:2015 on the electromagnetic compatibility of medical devices. Compliance with electromagnetic compatibility standards does not guarantee total immunity of the product. Some devices (mobile phones, pagers, etc.) may interrupt operation if used near medical devices.

9 - OPERATING AND STORAGE CONDITIONS

9.1 - OPERATING CONDITIONS

- temperature from +10 °C to +40 °C
- relative humidity from 0% to 95%

9.2 - STORAGE CONDITIONS

- temperature from -25 °C to +70 °C
- relative humidity from 10% a 95%
- atmospheric pressure from 690 hPa to 1060 hPa

10 - SPARE PARTS

Only original spare parts must be used.

11 - WARRANTY CONDITIONS

The manufacturer is liable for the safety, reliability and performance of the device, provided that it is used in accordance with the instructions and for the intended purpose only, and any repairs are carried out by the manufacturer or an approved service centre. The manufacturer declines all liability and the warranty will be invalidated if the user fails to follow the instructions and the important notes contained herein.

12 - TECHNICAL UPDATES

Air Liquide Medical Systems S.r.l. periodically reviews and modifies all its medical devices in order to improve their performance, safety and reliability. The instruction booklets are updated accordingly to include any new or changed features.

If the booklet accompanying the device gets damaged or mislaid, a replacement copy can be obtained from the manufacturer by quoting the data shown on the rating plate.

13 - COPYRIGHT

The information contained in this booklet must not be used for any purpose other than the reason for which it is provided.

This booklet is the property of Air Liquide Medical Systems S.r.l. and no part of it may be reproduced without the manufacturer's written permission. All rights reserved.

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